

ENDOACCESS

ENDO FOR LIFE



**INSTRUMENTS & ACCESSORIES
FOR FLEXIBLE ENDOSCOPES**



Disposable Instruments

Disposable ESD Instruments	Page 1 - 3
Disposable ERCP Instruments	Page 3 - 5
Disposable EMR Instruments	Page 5 - 8
Disposable Cold Snares	Page 8
Disposable Snare with Bag	Page 9
Disposable Hemoclips	Page 9
Disposable Grasping Forceps	Page 10 - 13
Disposable Retrieval Baskets	Page 14 - 15
Disposable Biopsy Forceps	Page 15 - 19
Disposable Hot Biopsy Forceps	Page 19
Disposable Spray Catheters	Page 20
Disposable Cytology Brushes	Page 21
Disposable Cleaning Brushes	Page 21 - 22
Disposable Accessories	Page 22 - 24

Reusable Instruments

Reusable Biopsy Forceps	Page 24 - 27
Reusable Foreign Body Removal	Page 27 - 31
Reusable Retrieval Baskets	Page 32 - 33
Reusable Electrical Electrodes	Page 33 - 34
Reusable Electrical Snares	Page 35 - 36
Reusable Scissors Type Grasping Forceps	Page 36 - 37
Reusable Cleaning Brushes	Page 37

ESD Instruments

Disposable Electrical Forceps (Rotatable, Hemo Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24HFR18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	7mm
EA-S24HFR23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7mm

mark, cut and stop bleeding by opening and closing the branches. Slow rotation facilitates operation
Olympus compatible HF connection

Disposable Electrical Forceps (Rotatable, Tooth/Hook Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24DTR18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	7mm
EA-S24DTR23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7mm

Easy cutting into the adhesion. Slow rotation facilitates fast access. Hook tip designed for lifting of tissue
Olympus compatible HF connection

Disposable Electrical Forceps (Rotatable, Dissection)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24DFR18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	7mm
EA-S24DFR23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7mm

Easy cutting into the adhesion. Slow rotation facilitates fast access.
Olympus compatible HF connection

Disposable Electrical Forceps (Rotatable, Dissection Scissors)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24DSR18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	7mm
EA-S24DSR23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7mm

Easy cutting into the adhesion. Slow rotation facilitates fast access
Olympus compatible HF connection

Disposable Ball Type



Item	Use	Tube Ø and Working Length	Working Channel Ø	Cutting Length
EA-S24BT18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	2mm
EA-S24BT23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	2mm

Multifunctional Design: marking, cutting, dissection, coagulation, spraying, etc

Ball Shape Type: Good operability, precise coagulation effect. Good traction during cutting and peeling

Spray function: For spraying into the submucosa to lift the mucosa or to show the bleeding area

Olympus compatible HF connection

Disposable Cylindrical Type (Tunnel Knife)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Cutting Length
EA-S24ST18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	1.5/2mm adjustable
EA-S24ST23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	1.5/2mm adjustable

Multifunctional Design: Marking, cutting, injection, lifting, dissection, coagulation, etc.

Needle shape tip: Sharp cutting tip extendable up to 2 mm

Submucosa injection function: High quality hydraulic solution pressure with needle point design; the knife can inject precisely into the submucosa and achieve a perfect mucosal lifting effect

Olympus compatible HF connection

Disposable Ball Type (CT Knife)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Cutting Length
EA-S24B3T18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	4mm
EA-S24B3T23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	4mm

Insulated Tip Design: Ensures safety when cutting and dissecting, excellent conductivity; gentle cutting and stable coagulation

Adjustable electrode length: Adjustable from 0 mm to 4 mm according to operation procedure and lesions

Spray function: For spraying into the submucosa to lift the mucosa or to show the bleeding area

Olympus compatible HF connection

Disposable Cylindrical Knife (Square Knife)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Cutting Length
EA-S24S2T18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	4mm
EA-S24S2T23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	4mm

Square Shape Tip Design: Marking and coagulation

The design of the tip allows to lift or hook the tissue to cut efficiently. Good support surface and better coagulation effect


Adjustable electrode length: Adjustable cutting tip from 0 to 4mm

Spray function: For spraying into the submucosa to lift the mucosa or to show the bleeding area

Olympus compatible HF connection

Disposable ERCP Instruments

Disposable Papillotome



Item	Tube Ø and Working Length	Working Channel Ø	Tip length A	Cutting length B
Double Lumen				
EA-S0720EPM2	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm	20mm
EA-S0725EPM2	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm	25mm
EA-S0730EPM2	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm	30mm
Triple Lumen				
EA-S0720EPM3	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm	20mm
EA-S0725EPM3	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm	25mm
EA-S0730EPM3	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm	30mm

braided wire, radiopaque markers ; preformed and tapered tip to Ø1.9mm

Disposable Sphincterotome



Item	Tube Ø and Working Length	Working Channel Ø	Needle Length
EA-S24NE212 (Double Lumen)	Ø2.4mm x 2100mm	≥Ø2.8mm	3-5mm (adjustable)
EA-S24NE213 (Triple Lumen)	Ø2.4mm x 2100mm	≥Ø2.8mm	3-5mm (adjustable)

adjustable needle length 3-5mm, monofil, radiopaque marker

Disposable ERCP Basket



Item	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24BS21-1 (Single Lumen)	Ø2.4mm x 2100mm	≥Ø2.8mm	20mm
EA-S24BS21-2 (Double Lumen)	Ø2.4mm x 2100mm	≥Ø2.8mm	20mm

For stone extraction, suitable for lithotripsy, hexagonal, retractable atraumatic tip, double lumen, basket length 45mm

Disposable ERCP Cytology Brush



Item	Tube Ø and Working Length	Working Channel Ø
EA-S24CB21-2 double lumen, brush length 10mm	Ø2.4mm x 2100mm	≥Ø2.8mm

Disposable Guide Wire



Item	Tube Ø and Working Length	Tip
EA-FUH25460	Ø0.025" x 4600mm	55 mm hydrophil/made in Germany
EA-FUH35460	Ø0.035" x 4600mm	55 mm hydrophil/made in Germany
EA-FUH35450	Ø0.035" x 4500mm	55 mm hydrophil

Disposable EMR Instruments

Disposable Injection Needle




- 1) stainless steel needle point, teflon inner & outer tube
- 2) EA-SxxxxNJxx xx

			-----needle length
			----- needle diameter
			----- tube diameter & working length

Item	Use	Tube Ø and Working Length	Working Channel Ø	Needle Length	Needle Ø
EA-S1810NJxxxx	Gastroscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	4/5mm	23/25G
EA-S1816NJxxxx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	4/5/6mm	23/25G
EA-S1825NJxxxx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	4/5/6mm	23/25G
EA-S2416NJxxxx	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	4/5/6mm	21/23/25G
EA-S2418NJxxxx	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	4/5/6mm	21/23/25G
EA-S2421NJxxxx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	4/5/6mm	21/23/25G
EA-S2423NJxxxx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	4/5/6mm	21/23/25G

with locking mechanism

Disposable Electrical Snare (Oval Type) multifil rotatable / non-rotatable



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
Rotatable				
EA-S18OS10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18OS12xx	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	10/15/20mm
EA-S18OS16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S18OS25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-S24OS16xx	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24OS18xx	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24OS21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24OS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

Non-rotatable				
EA-S18OSN10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18OSN12xx	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	10/15/20mm
EA-S18OSN16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S18OSN25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-S24OSN16xx	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24OSN18xx	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24OSN21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24OSN23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

Olympus compatible HF connection, infinitely adjustable

Disposable Electrical Snare (Oval Type) monofil rotatable / non-rotatable



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
Rotatable				
EA-S24OSM23xxR	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm
Non-rotatable				
EA-S24OSM23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

Olympus compatible HF connection, infinitely adjustable


Disposable Electrical Snare (Hexagonal Type) rotatable / non-rotatable



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
Rotatable				
EA-S18HS10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18HS12xx	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	10/15/20mm
EA-S18HS16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
Ea-S18HS25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-S24HS16xx	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24HS18xx	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24HS21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24HS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm
Non-rotatable				
EA-S18HSN10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18HSN12xx	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	10/15/20mm
EA-S18HSN16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S18HSN25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-S24HSN16xx	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24HSN18xx	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24HSN21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24HSN23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

Olympus compatible HF connection

Disposable Electrical Snare (Half-Moon Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18HM10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18HM12xx	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	10/15/20mm
EA-S18HM16xx	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	10/15/20mm
Ea-S18HM25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-S24HM16xx	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	20/25/30/mm
EA-S24HM18xx	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	20/25/30/mm
EA-S24HM21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	20/25/30/mm
EA-S24HM23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	20/25/30/mm

Asymmetric, Olympus compatible HF connection

Disposable Cold Snares

Disposable Grasping Forceps (Cold Snare Type, monofil)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24CS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15mm

Disposable Snare (Cold Snare Type with HF connection, monofil)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24OSC23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15mm

Olympus compatible HF connection

Disposable Snare with Bag

Disposable Snare with Bag



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18SN10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	20mm
EA-S18SN12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	20mm
EA-S18SN16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	20mm
EA-S24SN18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	20/30mm
EA-S24SN23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	20/30mm

Tear proof, perforated extraction bag for removal of foreign bodies and polyps

Disposable Hemoclip

Disposable Hemoclip



Item	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S24HC2309R	Ø2.4mm x 2300mm	≥Ø2.8mm	9mm
EA-S24HC2311R	Ø2.4mm x 2300mm	≥Ø2.8mm	11mm
EA-S24HC2314R	Ø2.4mm x 2300mm	≥Ø2.8mm	14mm
EA-S24HC2317R	Ø2.4mm x 2300mm	≥Ø2.8mm	17mm

Please ask for more available sizes

Preloaded, repositionable, 360° rotation

Disposable Grasping Forceps

Disposable Grasping Forceps (Rat Tooth Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18RG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5mm
EA-S18RG12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5mm
EA-S18RG16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	5mm
EA-S24RG16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6mm
EA-S24RG18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6mm
EA-S24RG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6mm
EA-S24RG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6mm

For the removal of flat and hard object (e.g. coins)

Disposable Grasping Forceps (Rat Tooth with Alligator Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18RAG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5mm
EA-S18RAG12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5mm
EA-S18RAG16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	5mm
EA-S24RAG16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6mm
EA-S24RAG18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6mm
EA-S24RAG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6mm
EA-S24RAG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6mm

Sharp teeth, long jaws for difficult retrieval of flat and hard objects

Disposable Grasping Forceps (Pelican Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18PG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	9.5mm
EA-S18PG12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	9.5mm
EA-S18PG16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	9.5mm
EA-S24PG16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	15mm
EA-S24PG18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	15mm
EA-S24PG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	15mm
EA-S24PG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	15mm

For the retrieval of big objects or stent removal

Disposable Grasping Forceps (Alligator Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18AG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	7mm
EA-S18AG12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	7mm
EA-S18AG16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	7mm
EA-S24AG16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	7.5mm
EA-S24AG18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	7.5mm
EA-S24AG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	7.5mm
EA-S24AG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7.5mm

For the removal of sharp-edged objects

Disposable Grasping Forceps (3-Prong Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S18TP07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-S18TP10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-S18TP12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm
EA-S18TP15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm
EA-S24TP16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-S24TP17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-S24TP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-S24TP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Round grippers, for the removal of big polyps and foreign bodies


Disposable Grasping Forceps (3-Prong Hook Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S18TP07H	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-S18TP10H	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-S18TP15H	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm

This item is also available on request in Ø 2.4. Longer lead times have to be considered.


Disposable Grasping Forceps (5-Prong Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S18FP07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-S18FP10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-S18FP15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm
EA-S24FP16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-S24FP17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-S24FP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-S24FP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Round grippers, for the removal of big polyps and foreign bodies

Disposable Grasping Forceps (5-Prong Hook Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S18FP07H	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-S18FP10H	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-S18FP15H	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm

This item is also available on request in Ø 2.4. Longer lead times have to be considered

Disposable Grasping Forceps (2-Prong Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S24DP16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-S24DP17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-S24DP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-S24DP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

For the removal of gastric balloons

Disposable Retrieval Baskets

Disposable Retrieval Basket (Fold Angular, stainless steel)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18BS07xx	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	10/15/20mm
EA-S18BS10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18BS16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S18BS25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-S24BS15xx	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24BS17xx	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24BS21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24BS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

For the removal of stones and round and slippery objects, multifilament wire, unsuitable for lithotripsy

Disposable Retrieval Basket (without Tip, Nitinol)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S10UB07xx	Urology	Ø1.0mm x 700mm	≥Ø1.2mm	10/15/mm
EA-S10UB10xx	Urology	Ø1.0mm x 1000mm	≥Ø1.2mm	10/15/mm
EA-S10UB12xx	Urology	Ø1.0mm x 1200mm	≥Ø1.2mm	10/15/mm
EA-S18BM07xx	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	10/15/20mm
EA-S18BM10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18BM16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S18BM25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-S24BM15xx	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24BM17xx	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24BM21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24BM23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

For the removal of stones and round and slippery objects, multifilament wire, unsuitable for lithotripsy

Disposable Spirality Basket (Nitinol)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S10U2B07xx	Urology	Ø1.0mm x 700mm	≥Ø1.2mm	10/15/mm
EA-S10U2B10xx	Urology	Ø1.0mm x 1000mm	≥Ø1.2mm	10/15/mm
EA-S10U2B12xx	Urology	Ø1.0mm x 1200mm	≥Ø1.2mm	10/15/mm
EA-S18B2M07xx	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	10/15/20mm
EA-S18B2M10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-S18B2M16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-S18B2M25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-S24B2M15xx	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24B2M17xx	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24B2M21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-S24B2M23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

For the removal of stones and round and slippery objects, multifilament wire, unsuitable for lithotripsy

Disposable Biopsy Forceps

Disposable Biopsy Forceps (Oval, fenestrated cups)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18C10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S18C12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S18C15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S18C16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	5.5mm
EA-S24C16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24C18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24C21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24C23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm
EA-S30C23	Colonoscopy	Ø3.0mm x 2300mm	≥Ø3.2mm	6.5mm

Oval, fenestrated cups ensure minimal damage of tissue. Allows deeper biopsies. Uncoated

Disposable Biopsy Forceps (Oval, coated, fenestrated cups)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S180C10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S180C12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S180C15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S180C16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	5.5mm
EA-S180C18	Gastroscopy	Ø1.8mm x 1800mm	≥Ø2.0mm	5.5mm
EA-S240C16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S240C18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S240C21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S240C23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm
EA-S300C23	Colonoscopy	Ø3.0mm x 2300mm	≥Ø3.2mm	6.5mm

Oval, fenestrated cups ensure minimal damage of tissue. Allows deeper biopsies. Coated

Disposable Biopsy Forceps (Oval Cups with Needle, fenestrated)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S180W10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S180W12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S180W15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S240W16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.0mm	6.5mm
EA-S240W18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S240W21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S240W23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Oval, fenestrated cups with needle. Improved fixation within tissue to avoid slipping. Minimal damage of tissue. Uncoated

Disposable Biopsy Forceps (Oval Cups with Needle, coated, fenestrated)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18CW10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S18CW12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S18CW15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S24CW16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24CW18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24CW21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24CW23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Oval, fenestrated cups with needle. Improved fixation within tissue to avoid slipping. Minimal damage of tissue. Coated

Disposable Biopsy Forceps (Alligator Cups)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18A10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S18A12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S18A15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S24A16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24A18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24A21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24A23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Serrated alligator cups for deep biopsy of hard and difficult to access tissue. Uncoated

Disposable Biopsy Forceps (Alligator Cups, coated)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18AC10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S18AC12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S18AC15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S24AC16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24AC18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24AC21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24AC23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Serrated alligator cups for deep biopsy of hard and difficult to access tissue. Coated

Disposable Biopsy Forceps (Alligator Cups with Needle, fenestrated)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18AW10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S18AW12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S18AW15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S24AW16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24AW18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24AW21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24AW23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Serrated alligator cups for deep biopsy of hard and difficult to access tissue. Improved fixation within tissue. Uncoated

Disposable Biopsy Forceps (Alligator cups with Needle, coated)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18WA10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-S18WA12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5.5mm
EA-S18WA15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-S24WA16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24WA18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24WA21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24WA23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Serrated alligator cups for deep biopsy of hard and difficult to access tissue. Improved fixation within tissue.
Coated

Disposable Hot Biopsy Forceps

Disposable Hot Biopsy Forceps



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-S18H10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5mm
EA-S18H12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	5mm
EA-S18H15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5mm
EA-S24H16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	6.5mm
EA-S24H18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	6.5mm
EA-S24H21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-S24H23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

For the ablation of small lesions with simultaneous hemostasis or ablation of diminutive polyps. Olympus compatible HF connection

Disposable Spray Catheters

Disposable Standard Spray Catheter



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S18SC10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-S18SC12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm
EA-S18SC16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm
EA-S24SC16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-S24SC18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm
EA-S24SC21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-S24SC23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Even spray for clear images. For dye, saline and other solutions

Disposable Fog Spray Catheters



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S18FC10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-S18FC12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm
EA-S18FC16	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm
EA-S24FC16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-S24FC18	Gastroscopy	Ø2.4mm x 1800mm	≥Ø2.8mm
EA-S24FC21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-S24FC23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Extensive spray. For dye, saline and other solutions

Disposable Cytology Brushes

Disposable Cytology Brushes



Item	Use	Tube Ø and Working Length	Working Channel Ø	Ø Brushes and Length
EA-S18CB10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	3mm/7mm
EA-S18CB12	Bronchoscopy	Ø1.8mm x 1200mm	≥Ø2.0mm	3mm/7mm
EA-S18CB15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	3mm/7mm
EA-S24CB16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm	3mm/7mm
EA-S24CB18	Colonoscopy	Ø2.4mm x 1800mm	≥Ø2.8mm	3mm/7mm
EA-S24CB23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	3mm/7mm

Excellent brush hardness for improved sample collection. Thinner bristles for significantly improved cell collection capability.

Disposable Cleaning Brushes


Disposable Cleaning Brush (Double ended)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-S10DB10	Urology	Ø1.0mm x 1000mm	≥Ø1.2mm
EA-S10DB12	Urology	Ø1.0mm x 1200mm	≥Ø1.2mm
EA-S10DB22	Urology	Ø1.0mm x 2200mm	≥Ø1.2mm
EA-S18DB18	Bronchoscopy	Ø1.8mm x 1800mm	≥Ø2.0mm
EA-S24DB16	Colonoscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-S24DB22	Colonoscopy	Ø2.4mm x 2200mm	≥Ø2.8mm

For cleaning the working channel, symmetrical

Disposable Cleaning Brush (Double ended, asymmetrical)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Ø Brushes
EA-S18DB18A34	Bronchoscopy	Ø1.8mm x 1800mm	≥Ø2.0mm	3mm/4mm
EA-S18DB23A35	Bronchoscopy	Ø1.8mm x 2300mm	≥Ø2.0mm	3mm/5mm
EA-S24DB22A510	Colonoscopy	Ø2.4mm x 2200mm	≥Ø2.8mm	5mm/10mm

For cleaning the working channel, asymmetrical

Disposable Port Hole Cleaning Brush (Double ended, asymmetrical)

Valve cleaning



Item	Ø Brushes
EA-S1001PCD	5/10mm

Disposable Port Hole Cleaning Brush

Valve cleaning



Item	Ø Brush
EA-S1001PC	10mm

Disposable Accessories

Disposable Colonoscopy Pants



Item	Size
EA-SCP-L	L (blue)
EA-SCP-XL	XL (green)

Polyp Trap (4 chambers)



Item	Chambers
EA-SPT Removable chambers, transparent, rotatable	4
EA-SPT-S Transparent	1

Disposable Helico-Check



Item
EA-Prontrodry
Disposable helico check for the detection of H. pylori / results within 5 minutes

Disposable Bite Block



Item	Strap
EA-SBBW	Yes
For adults	

Disposable Bite Block (edentulous patients, children)



Item EA-SBBK-N	Strap Yes
--------------------------	---------------------

Disposable Bite Block (children)



Item EA-SBBK	Strap Yes
------------------------	---------------------

Disposable Valve Cap



Item EA-SVC	Compatibility Olympus
-----------------------	---------------------------------

Closable lid

Disposable Valve Cap Set



Item EA-SVC-S	Compatibility Olympus
-------------------------	---------------------------------

Valve cap, air/water valve; suction valve

Reusable Biopsy Forceps

Reusable Biopsy Forceps (Standard Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18S07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5mm
EA-R18S10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5mm
EA-R18S15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5mm
EA-R18S25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	5mm
EA-R24S15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	5.5mm
EA-R24S17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	5.5mm
EA-R24S21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	5.5mm
EA-R24S23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	5.5mm

Standard, round fenestrated cups for minimal tissue damage, autoclavable

Reusable Biopsy Forceps (Oval Fenestrated Cups Type)



Use	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18C07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5.5mm
EA-R18C10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18C15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R18C25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	5.5mm
EA-R24C15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	6.5mm
EA-R24C17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	6.5mm
EA-R24C21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-R24C23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Oval, fenestrated cups for deeper biopsy, autoclavable

Reusable Biopsy Forceps (Oval Fenestrated Cups with Needle)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R180W07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5mm
EA-R180W10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5mm
EA-R180W15	Gastrosocopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5mm
EA-R180W25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	5mm
EA-R240W15	Gastrosocopy	Ø2.4mm x 1500mm	≥Ø2.8mm	7.5mm
EA-R240W17	Gastrosocopy	Ø2.4mm x 1700mm	≥Ø2.8mm	7.5mm
EA-R240W21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	7.5mm
EA-R240W23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7.5mm

Oval, fenestrated cups with needle. Improved fixation within tissue to avoid slipping.
Minimal damage of tissue, autoclavable

Reusable Biopsy Forceps (Alligator Jaws Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18A07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5.5mm
EA-R18A10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18A15	Gastrosocopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R18A25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	5.5mm
EA-R24A15	Gastrosocopy	Ø2.4mm x 1500mm	≥Ø2.8mm	6.5mm
EA-R24A17	Gastrosocopy	Ø2.4mm x 1700mm	≥Ø2.8mm	6.5mm
EA-R24A21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-R24A23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Serrated alligator cups for deep biopsy of hard and difficult to access tissue, autoclavable

Reusable Biopsy Forceps (Alligator, Single Moving Cup Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18OP07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5.5mm
EA-R18OP10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18OP15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R18OP25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	5.5mm
EA-R24OP15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	6.5mm
EA-R24OP17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	6.5mm
EA-R24OP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-R24OP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Single tilting cup (90°) for tangential sampling, autoclavable

Reusable Biopsy Forceps (Alligator Jaws with Needle)

Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18AW07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5.5mm
EA-R18AW10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18AW15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R18AW25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	5.5mm
EA-R24AW15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	6.5mm
EA-R24AW17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	6.5mm
EA-R24AW21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	6.5mm
EA-R24AW23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	6.5mm

Serrated alligator cups for deep biopsy of hard and difficult to access tissue. Improved fixation within tissue.
Uncoated, autoclavable

Reusable Foreign Body Grasping Forceps

Reusable Grasping Forceps (Pelican Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18PG07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	11.5mm
EA-R18PG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	11.5mm
EA-R18PG15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	11.5mm
EA-R18PG25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	11.5mm
EA-R24PG15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	15mm
EA-R24PG17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	15mm
EA-R24PG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	15mm
EA-R24PG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	15mm

For the retrieval of big objects or stent removal, autoclavable

Reusable Grasping Forceps (Alligator Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18AG07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	7mm
EA-R18AG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	7mm
EA-R18AG15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	7mm
EA-R18AG25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	7mm
EA-R24AG15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	8mm
EA-R24AG17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	8mm
EA-R24AG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	8mm
EA-R24AG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	8mm

For the removal of sharp-edged objects, autoclavable

Reusable Grasping Forceps (Rat Tooth Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18RG07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5.5mm
EA-R18RG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18RG15	Gastrosocopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R18RG25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	5.5mm
EA-R24RG15	Gastrosocopy	Ø2.4mm x 1500mm	≥Ø2.8mm	5.5mm
EA-R24RG17	Gastrosocopy	Ø2.4mm x 1700mm	≥Ø2.8mm	5.5mm
EA-R24RG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	5.5mm
EA-R24RG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	5.5mm

For the removal of flat and hard object (e.g. coins), autoclavable

Reusable Grasping Forceps (Rat Tooth Type with Alligator Jaws Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18RAG07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5.5mm
EA-R18RAG10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18RAG15	Gastrosocopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R18RAG25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	5.5mm
EA-R24RAG15	Gastrosocopy	Ø2.4mm x 1500mm	≥Ø2.8mm	5.5mm
EA-R24RAG17	Gastrosocopy	Ø2.4mm x 1700mm	≥Ø2.8mm	5.5mm
EA-R24RAG21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	5.5mm
EA-R24RAG23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	5.5mm

Sharp teeth, long jaws for difficult retrieval of flat and hard objects, autoclavable


Reusable Grasping Forceps (3-Prong Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R18TP07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-R18TP10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-R18TP15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm
EA-R24TP16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-R24TP17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-R24TP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-R24TP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Round grippers, for the removal of big polyps and foreign bodies, autoclavable


Reusable Grasping Forceps (3-Prong Hook Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R18TP07H	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-R18TP10H	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-R18TP15H	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm

autoclavable

Reusable Grasping Forceps (5-Prong Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R18FP07	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-R18FP10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-R18FP15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm
EA-R24FP16	Gastroscopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-R24FP17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-R24FP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-R24FP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm


Round grippers, for the removal of big polyps and foreign bodies, autoclavable

Reusable Grasping Forceps (5-Prong HookType)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R18FP07H	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm
EA-R18FP10H	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-R18FP15H	Gastrosocopy	Ø1.8mm x 1500mm	≥Ø2.0mm
EA-R24FP23H autoclavable	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Reusable Grasping Forceps (2-Prong Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R24DP16	Gastrosocopy	Ø2.4mm x 1600mm	≥Ø2.8mm
EA-R24DP17	Gastrosocopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-R24DP21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-R24DP23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

For the removal of gastric balloons, autoclavable

Reusable Grasping Forceps (Rubber Tips Type)




Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18RT10	Bronchoscopy	Ø 1.8mm x 1000mm	≥Ø2.0mm	5.5mm
EA-R18RT15	Gastrosocopy	Ø 1.8mm x 1500mm	≥Ø2.0mm	5.5mm
EA-R24RT15	Gastrosocopy	Ø 2.4mm x 1500mm	≥Ø2.8mm	6mm
EA-R24RT17	Gastrosocopy	Ø 2.4mm x 1700mm	≥Ø2.8mm	6mm
EA-R24RT21	Colonoscopy	Ø 2.4mm x 2100mm	≥Ø2.8mm	6mm
EA-R24RT23	Colonoscopy	Ø 2.4mm x 2300mm	≥Ø2.8mm	6mm

For grasping slippery, flat and sharpedged objects, e.g. needles, wire, autoclavable

Reusable Retrieval Baskets

Reusable Fold Angular Basket (Stainless Steel)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18BS07xx	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	10/15/20mm
EA-R18BS10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-R18BS16xx	Gastrosocopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-R18BS25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-R24BS15xx	Gastrosocopy	Ø2.4mm x 1500mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24BS17xx	Gastrosocopy	Ø2.4mm x 1700mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24BS21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24BS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

For the removal of stones and round and slippery objects, multifilament wire, unsuitable for lithotripsy, autoclavable


Reusable Basket without Tip (Nitinol)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18BM07xx	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	10/15/20mm
EA-R18BM10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-R18BM16xx	Gastrosocopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-R18BM25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-R24BM15xx	Gastrosocopy	Ø2.4mm x 1500mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24BM17xx	Gastrosocopy	Ø2.4mm x 1700mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24BM21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24BM23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

For the removal of stones and round and slippery objects, multifilament wire, unsuitable for lithotripsy, autoclavable

Reusable Spirality Basket (Nitinol)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18SB07xx	Choledochoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	10/15/20mm
EA-R18SB10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-R18SB16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-R18SB25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-R24SB15xx	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24SB17xx	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24SB21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24SB23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

For the removal of stones and round and slippery objects, multifilament wire, unsuitable for lithotripsy, autoclavable

Reusable Electrical Electrodes

Reusable Electrical Forceps (Coagulation Electrode)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R18CE10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-R18CE15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm
EA-R18CE25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm
EA-R24CE15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm
EA-R24CE17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-R24CE21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-R24CE23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Coagulation of blood vessels, autoclavable

Reusable Electrical Forceps (Hot Biopsy)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18H10	Bronchoscopy	Ø1.8mm x 700mm	≥Ø2.0mm	5mm
EA-R18H15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm	5mm
EA-R18H25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	5mm
EA-R24H15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.0mm	7mm
EA-R24H17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm	7mm
EA-R24H21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	7mm
EA-R24H23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	7mm

For the ablation of small lesions with simultaneous hemostasis or ablation of diminutive polyps, autoclavable

Reusable Electrical Forceps (Needle Electrode)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R18NE10	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm
EA-R18NE15	Gastroscopy	Ø1.8mm x 1500mm	≥Ø2.0mm
EA-R18NE25	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm
EA-R24NE15	Gastroscopy	Ø2.4mm x 1500mm	≥Ø2.8mm
EA-R24NE17	Gastroscopy	Ø2.4mm x 1700mm	≥Ø2.8mm
EA-R24NE21	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm
EA-R24NE23	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm

Cutting tissue with high frequency; needle length 4mm, autoclavable

Reusable Electrical Snares


Reusable Electrical Snare (Oval Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18OS10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-R18OS16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-R18OS25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-R24OS16xx	Gastroscopy	Ø 2.4mm x 1600mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24OS21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24OS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

Olympus compatible HF connection, infinitely adjustable, autoclavable

Reusable Electrical Snare (Hexagonal Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18HS10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-R18HS16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-R18HS25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-R24HS16xx	Gastroscopy	Ø 2.4mm x 1600mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24HS21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	10/15/20/25/30/35mm
EA-R24HS23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	10/15/20/25/30/35mm

Olympus compatible HF connection, autoclavable

Reusable Electrical Snare (Half-Moon Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R18HM10xx	Bronchoscopy	Ø1.8mm x 1000mm	≥Ø2.0mm	10/15/20mm
EA-R18HM16xx	Gastroscopy	Ø1.8mm x 1600mm	≥Ø2.0mm	10/15/20mm
EA-R18HM25xx	Enteroscopy	Ø1.8mm x 2500mm	≥Ø2.0mm	10/15/20mm
EA-R24HM16xx	Gastroscopy	Ø 2.4mm x 1600mm	≥Ø2.8mm	20/25/30mm
EA-R24HM21xx	Colonoscopy	Ø2.4mm x 2100mm	≥Ø2.8mm	20/25/30mm
EA-R24HM23xx	Colonoscopy	Ø2.4mm x 2300mm	≥Ø2.8mm	20/25/30mm

Asymmetric, Olympus compatible HF connection, autoclavable

Reusable Scissors Type Grasping Forceps

Reusable Grasping Forceps (Scissors Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø	Opening
EA-R24ST15	Gastroscopy	Ø 2.4mm x 1500mm	≥Ø2.8mm	6mm
EA-R24ST17	Gastroscopy	Ø 2.4mm x 1700mm	≥Ø2.8mm	6mm
EA-R24ST21	Colonoscopy	Ø 2.4mm x 2100mm	≥Ø2.8mm	6mm
EA-R24ST23	Colonoscopy	Ø 2.4mm x 2300mm	≥Ø2.8mm	6mm

For cutting tissue, autoclavable

Reusable Grasping Forceps (Half-Scissors Type)



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R24HT15	Gastroscopy	Ø 2.4mm x 1500mm	≥Ø2.8mm
EA-R24HT17	Gastroscopy	Ø 2.4mm x 1700mm	≥Ø2.8mm
EA-R24HT21	Colonoscopy	Ø 2.4mm x 2100mm	≥Ø2.8mm
EA-R24HT23	Colonoscopy	Ø 2.4mm x 2300mm	≥Ø2.8mm

For cutting protruding tissue after ligature, autoclavable

Reusable Cleaning Brushes

Reusable Cleaning Brush



Item	Use	Tube Ø and Working Length	Working Channel Ø
EA-R18D16	Gastroscopy	Ø 1.8mm x 1600mm	≥Ø2.0mm
EA-R24D16	Colonoscopy	Ø 2.4mm x 1600mm	≥Ø2.8mm
EA-R24D18	Colonoscopy	Ø 2.4mm x 1800mm	≥Ø2.8mm
EA-R24D22	Colonoscopy	Ø 2.4mm x 2200mm	≥Ø2.8mm

For cleaning the working channel,

Reusable Porthole Cleaning Brush




Item
EA-R10PC


Valve cleaning

The logo consists of a yellow circle partially overlapping a black stylized shape that resembles a lowercase 'g' or a curved line.

ENDOACCESS
ENDO FOR LIFE

A close-up photograph of a silver-colored metal surgical instrument with two serrated jaws, positioned in the upper left quadrant.

**Endoaccess GmbH
Feldriethe 1
30826 Garbsen
Germany**

A close-up photograph of a silver-colored metal surgical instrument, similar to the one in the top left, positioned in the bottom right corner.

**+49 5131-44 22 61 0
info@endoaccess.de
www.endoaccess.de**