



Accessories for Flexible Endoscopy

Dear Sir or Madam,

Thank you for your interest in our products.

Endoaccess GmbH with its headquarters in Garbsen, Germany, is a producer and supplier of accessories for flexible endoscopes. We have been supplying our instruments to physicians of human medicine for many years now and are happy to now being able to deliver to veterinary practioners as well.

Our high quality instruments will enable you to treat your patients reliably at low costs.

All sizes displayed are standard sizes.

We also offer custom-made production. Please inquire for possibilites and prices.

All reusable items are also available as single-use products.

Foreign Body Removal: Grasping Forceps	1 - 5
Foreign Body Removal: Snares	6 - 8
Foreign Body Removal: Retrieval Baskets	9 - 11
Electrosurgical Instruments (HF)	12- 13
Biopsy Forceps	14- 16
Accessories	17- 19

1

Reusable Grasping Forceps


Rat Tooth



For the removal of flat or hard objects
opening 5.5mm (Ø 1.8mm & Ø 2.4mm)

- EA-R18RG10 Ø 1.8mm x 1000mm
- EA-R18RG12 Ø 1.8mm x 1200 mm
- EA-R18RG15 Ø 1.8mm x 1500mm

- EA-R24RG17 Ø 2.4mm x 1700mm
- EA-R24RG18 Ø 2.4mm x 1800mm
- EA-R24RG21 Ø 2.4mm x 2100mm
- EA-R24RG23 Ø 2.4mm x 2300mm




Alligator



For the removal of sharp-edged objects
opening 7mm (Ø 1.8mm)
8mm (Ø 2.4mm)

- EA-R18AG10 Ø 1.8mm x 1000mm
- EA-R18AG12 Ø 1.8mm x 1200 mm
- EA-R18AG15 Ø 1.8mm x 1500mm

- EA-R24AG17 Ø 2.4mm x 1700mm
- EA-R24AG18 Ø 2.4mm x 1800mm
- EA-R24AG21 Ø 2.4mm x 2100mm
- EA-R24AG23 Ø 2.4mm x 2300mm



2

Reusable Grasping Forceps

**Rat Tooth with
Alligator Jaws**



**Sharp teeth and long jaws for difficult
removal of flat and hard objects**

opening 5.5mm (Ø 1.8mm & Ø 2.4mm)

EA-R18RAG10 Ø 1.8mm x 1000mm
EA-R18RAG12 Ø 1.8mm x 1200 mm
EA-R18RAG15 Ø 1.8mm x 1500mm

EA-R24RAG17 Ø 2.4mm x 1700mm
EA-R24RAG18 Ø 2.4mm x 1800mm
EA-R24RAG21 Ø 2.4mm x 2100mm
EA-R24RAG23 Ø 2.4mm x 2300mm




Pelican



**For the removal of big objects
opening 11.5mm (Ø 1.8mm)
15mm (Ø 2.4mm)**


EA-R18PG10 Ø 1.8mm x 1000mm
EA-R18PG12 Ø 1.8mm x 1200 mm
EA-R18PG15 Ø 1.8mm x 1500mm

EA-R24PG17 Ø 2.4mm x 1700mm
EA-R24PG18 Ø 2.4mm x 1800mm
EA-R24PG21 Ø 2.4mm x 2100mm
EA-R24PG23 Ø 2.4mm x 2300mm



3



Reusable Grasping Forceps



For the removal of big polyps and foreign bodies

EA-R18TP10 Ø 1.8mm x 1000mm
EA-R18TP12 Ø 1.8mm x 1200 mm
EA-R18TP15 Ø 1.8mm x 1500mm


EA-R24TP17 Ø 2.4mm x 1700mm
EA-R24TP18 Ø 2.4mm x 1800mm
EA-R24TP21 Ø 2.4mm x 2100mm
EA-R24TP23 Ø 2.4mm x 2300mm



For the removal of big polyps and foreign bodies
Suitable for soft, brittle tissue

EA-R18FP10 Ø 1.8mm x 1000mm
EA-R18FP12 Ø 1.8mm x 1200 mm
EA-R18FP15 Ø 1.8mm x 1500mm

EA-R24FP17 Ø 2.4mm x 1700mm
EA-R24FP18 Ø 2.4mm x 1800mm
EA-R24FP21 Ø 2.4mm x 2100mm
EA-R24FP23 Ø 2.4mm x 2300mm



4

Reusable Grasping Forceps Scissors Type




Rubber Tips

For grasping slippery, flat and sharp-edged objects, e.g. needles, wire

opening 5.5mm (Ø 1.8mm)
6mm (Ø 2.4mm)

- EA-R18RT10 Ø 1.8mm x 1000mm
- EA-R18RT12 Ø 1.8mm x 1200mm
- EA-R18RT15 Ø 1.8mm x 1500mm

- EA-R24RT15 Ø 2.4mm x 1500mm
- EA-R24RT17 Ø 2.4mm x 1700mm
- EA-R24RT21 Ø 2.4mm x 2100mm
- EA-R24RT23 Ø 2.4mm x 2300mm




Scissors

For cutting tissue

opening 5.5mm (Ø 1.8mm)
6mm (Ø 2.4mm)

- EA-R18ST10 Ø 1.8mm x 1000mm
- EA-R18ST12 Ø 1.8mm x 1200mm
- EA-R18ST15 Ø 1.8mm x 1500mm

- EA-R24ST15 Ø 2.4mm x 1500mm
- EA-R24ST17 Ø 2.4mm x 1700mm
- EA-R24ST21 Ø 2.4mm x 2100mm
- EA-R24ST23 Ø 2.4mm x 2300mm



5 Reusable Grasping Forceps Scissors Type



Half-Scissors Type

**For cutting protruding tissue after
ligature**

- | | |
|-------------------|-------------------------|
| EA-R18HT10 | Ø 1.8mm x 1000mm |
| EA-R18HT12 | Ø 1.8mm x 1200mm |
| EA-R18HT15 | Ø 1.8mm x 1500mm |
| EA-R24HT15 | Ø 2.4mm x 1500mm |
| EA-R24HT17 | Ø 2.4mm x 1700mm |
| EA-R24HT21 | Ø 2.4mm x 2100mm |
| EA-R24HT23 | Ø 2.4mm x 2300mm |



6 Single-Use Snare with Bag



Easy removal of hard to grasp objects, e.g. marbles
Perforated bag enables drainage of fluids

Ø 1.8mm snare: opening 20mm

Ø 2.4mm snare: opening 30mm

EA-S18SN12 **Ø 1.8mm x 1200mm**
EA-S18SN16 **Ø 1.8mm x 1600mm**

EA-S24SN18 **Ø 2.4mm x 1800mm**
EA-S24SN23 **Ø 2.4mm x 2300mm**




7

Reusable Snares with HF connection

**Snares for the removal of foreign bodies and cutting polyps.
Optional HF connection.
HF: high degree of coagulation and precise cutting**

Set of snares contains 2 snares and 1 handle



Snare, oval, multifilament wire, HF connection

- EA-R18OS10 Ø 1.8mm x 1000mm**
- EA-R18OS15 Ø 1.8mm x 1500mm**
- EA-R18OS25 Ø 1.8mm x 2500mm**

**each available with these openings:
10/15/20mm**


- EA-R24OS15 Ø 2.4mm x 1500mm**
- EA-R24OS17 Ø 2.4mm x 1700mm**
- EA-R24OS21 Ø 2.4mm x 2100mm**
- EA-R24OS23 Ø 2.4mm x 2300mm**

**each available with these openings:
10/15/20/25/30/35mm**



8

Reusable Snare with HF connection



Snare, asymmetrical, multifilament wire,
HF connection

- EA-R18HM10 Ø 1.8mm x 1000mm
- EA-R18HM15 Ø 1.8mm x 1500mm
- EA-R18HM25 Ø 1.8mm x 2500mm

each available with these openings:
10/15/20mm

- EA-R24HM15 Ø 2.4mm x 1500mm
- EA-R24HM17 Ø 2.4mm x 1700mm
- EA-R24HM21 Ø 2.4mm x 2100mm
- EA-R24HM23 Ø 2.4mm x 2300mm

each available with these openings:
10/15/20/25/30/35mm



Snare, hexagonal, multifilament wire,
HF connection

- EA-R18HS10 Ø 1.8mm x 1000mm
- EA-R18HS15 Ø 1.8mm x 1500mm
- EA-R18HS25 Ø 1.8mm x 2500mm

each available with these openings:
10/15/20mm

- EA-R24HS15 Ø 2.4mm x 1500mm
- EA-R24HS17 Ø 2.4mm x 1700mm
- EA-R24HS21 Ø 2.4mm x 2100mm
- EA-R24HS23 Ø 2.4mm x 2300mm


each available with these openings:
10/15/20/25/30/35mm



9

Reusable Retrieval Baskets

spirality, Nitinol



Nitinol / with tip
Easy removal of foreign bodies

EA-R18SB0715
opening: 15mm

Ø 1.8mm x 700mm

EA-R18SB1020
EA-R18SB1920
opening: 20mm

Ø 1.8mm x 1000mm
Ø 1.8mm x 1900mm



EA-R24SB1720
opening: 20mm

Ø 2.4mm x 1700mm

EA-R24SB1925
EA-R24SB2325
opening: 25mm

Ø 2.4mm x 1900mm
Ø 2.4mm x 2300mm



EA-R24SB1530
EA-R24SB1930
EA-R24SB2130
EA-R24SB2330
opening: 30mm

Ø 2.4mm x 1500mm
Ø 2.4mm x 1900mm
Ø 2.4mm x 2100mm
Ø 2.4mm x 2300mm

EA-R24SB2335
opening: 35mm

Ø 2.4mm x 2300mm

10 Reusable Retrieval Baskets

foldable, Nitinol



Easy retrieval of round and smooth or slippery objects

EA-R18BM0715 **Ø 1.8mm x 700mm**
opening: 15mm

EA-R18BM1020 **Ø 1.8mm x 1000mm**
EA-R18BM1920 **Ø 1.8mm x 1900mm**
opening: 20mm

EA-R24BM1530 **Ø 2.4mm x 1500mm**
opening: 30mm


EA-R24BM1720 **Ø 2.4mm x 1700mm**
opening: 20mm

EA-R24BM2130 **Ø 2.4mm x 2100mm**
EA-R24BM2330 **Ø 2.4mm x 2300mm**
opening: 30mm



11 Reusable Retrieval Baskets

foldable, stainless steel



With tip / easy retrieval of round and smooth or slippery objects

EA-R18BS0715
opening: 15mm

Ø 1.8mm x 700mm

EA-R18BS1020
EA-R18BM1920
opening: 20mm

Ø 1.8mm x 1000mm
Ø 1.8mm x 1900mm



EA-R24BS1530
opening: 30mm

Ø 2.4mm x 1500mm

EA-R24BS1720
opening: 20mm

Ø 2.4mm x 1700mm



EA-R24BS2130
EA-R24BS2330
opening: 30mm

Ø 2.4mm x 2100mm
Ø 2.4mm x 2300mm

12 Reusable Elektrosurgical Instruments


Electrical Needle Electrode



Cutting tissue with high frequency
Needle length 4mm

EA-R18NE10 Ø 1.8mm x 1000mm
EA-R18NE15 Ø 1.8mm x 1500mm
EA-R18NE25 Ø 1.8mm x 2500mm

EA-R24NE15 Ø 2.4mm x 1500mm
EA-R24NE17 Ø 2.4mm x 1700mm
EA-R24NE21 Ø 2.4mm x 2100mm
EA-R24NE23 Ø 2.4mm x 2300mm



Hot Biopsy Forceps



**For the ablation of small lesions with
simultaneous hemostasis or ablation of
dimunitive polyps**
opening 5mm (Ø 1.8mm)
7mm (Ø 2.4mm)

EA-R18H10 Ø 1.8mm x 1000mm
EA-R18H15 Ø 1.8mm x 1500mm
EA-R18H25 Ø 1.8mm x 2500mm

EA-R24H15 Ø 2.4mm x 1500mm
EA-R24H17 Ø 2.4mm x 1700mm
EA-R24H21 Ø 2.4mm x 2100mm
EA-R24H23 Ø 2.4mm x 2300mm



13 Reusable Elektrosurgical Instruments

Electrical Coagulation Electrode



Coagulation of blood vessels

EA-R18CE10 Ø 1.8mm x 1000mm
EA-R18CE15 Ø 1.8mm x 1500mm
EA-R18CE25 Ø 1.8mm x 2500mm

EA-R24CE15 Ø 2.4mm x 1500mm
EA-R24CE17 Ø 2.4mm x 1700mm
EA-R24CE21 Ø 2.4mm x 2100mm
EA-R24CE23 Ø 2.4mm x 2300mm



Includes set of exchangeable tips

14 Reusable Biopsy Forceps

Standard




Round, fenestrated cups ensure minimal damage of vessels

**opening 5mm (Ø 1.8mm)
5.5mm (Ø 2.4mm)**

EA-R18S10 Ø 1.8mm x 1000mm
EA-R18S12 Ø 1.8mm x 1200mm
EA-R18S15 Ø 1.8mm x 1500mm

EA-R24S16 Ø 2.4mm x 1600mm
EA-R24S17 Ø 2.4mm x 1700mm
EA-R24S18 Ø 2.4mm x 1800mm
EA-R24S21 Ø 2.4mm x 2100mm
EA-R24S23 Ø 2.4mm x 2300mm



Oval




Oval, fenestrated cups for deeper biopsy

**opening 5.5mm (Ø 1.8mm)
6.5mm (Ø 2.4mm)**

EA-R18C10 Ø 1.8mm x 1000mm
EA-R18C12 Ø 1.8mm x 1200mm
EA-R18C15 Ø 1.8mm x 1500mm

EA-R24C16 Ø 2.4mm x 1600mm
EA-R24C17 Ø 2.4mm x 1700mm
EA-R24C18 Ø 2.4mm x 1800mm
EA-R24C21 Ø 2.4mm x 2100mm
EA-R24C23 Ø 2.4mm x 2300mm



15 Reusable Biopsy Forceps

Oval with
Needle



Oval, fenestrated cups with needle. Improved fixation within tissue to avoid slipping. Minimal damage of tissue

**opening 5mm (Ø 1.8mm)
7.5mm (Ø 2.4mm)**

EA-R18OW10 Ø 1.8mm x 1000mm
EA-R18OW12 Ø 1.8mm x 1200mm
EA-R18OW15 Ø 1.8mm x 1500mm

EA-R24OW16 Ø 2.4mm x 1600mm
EA-R24OW17 Ø 2.4mm x 1700mm
EA-R24OW18 Ø 2.4mm x 1800mm
EA-R24OW21 Ø 2.4mm x 2100mm
EA-R24OW23 Ø 2.4mm x 2300mm



Alligator



Standard alligator cups for taking tissues samples or removal of foreign bodies

**opening 5.5mm (Ø 1.8mm)
6.5mm (Ø 2.4mm)**

EA-R18A10 Ø 1.8mm x 1000mm
EA-R18A12 Ø 1.8mm x 1200mm
EA-R18A15 Ø 1.8mm x 1500mm

EA-R24A16 Ø 2.4mm x 1600mm
EA-R24A17 Ø 2.4mm x 1700mm
EA-R24A18 Ø 2.4mm x 1800mm
EA-R24A21 Ø 2.4mm x 2100mm
EA-R24A23 Ø 2.4mm x 2300mm



16 Reusable Biopsy Forceps


Alligator with Needle



Alligator cups with needle for improved fixation within tissue

**opening 5.5mm (Ø 1.8mm)
6.5mm (Ø 2.4mm)**

- | | |
|-------------------|-------------------------|
| EA-R18AW10 | Ø 1.8mm x 1000mm |
| EA-R18AW12 | Ø 1.8mm x 1200mm |
| EA-R18AW15 | Ø 1.8mm x 1500mm |
| EA-R24AW16 | Ø 2.4mm x 1600mm |
| EA-R24AW17 | Ø 2.4mm x 1700mm |
| EA-R24AW18 | Ø 2.4mm x 1800mm |
| EA-R24AW21 | Ø 2.4mm x 2100mm |
| EA-R24AW23 | Ø 2.4mm x 2300mm |




Alligator, single tilting 90°



Single tilting cup for tangential entry


- | | |
|-------------------|-------------------------|
| EA-R18OP10 | Ø 1.8mm x 1000mm |
| EA-R18OP12 | Ø 1.8mm x 1200mm |
| EA-R18OP15 | Ø 1.8mm x 1500mm |
| EA-R24OP16 | Ø 2.4mm x 1600mm |
| EA-R24OP17 | Ø 2.4mm x 1700mm |
| EA-R24OP18 | Ø 2.4mm x 1800mm |
| EA-R24OP21 | Ø 2.4mm x 2100mm |
| EA-R24OP23 | Ø 2.4mm x 2300mm |



17

Accessories

Cleaning Brushes



**Double-ended channel-cleaning brush
single-use**

EA-S10DB12	Ø 1.0mm x 1200mm
EA-S10DB22	Ø 1.0mm x 2200mm
EA-S18DB18	Ø 1.8mm x 1800mm
EA-S24DB22	Ø 2.4mm x 2200mm


unit = 10 pcs



**Single-ended channel-cleaning brush
reusable**

EA-R18D16	Ø 1.8mm x 1600mm
EA-R24D16	Ø 2.4mm x 1600mm
EA-R24D18	Ø 2.4mm x 1800mm
EA-R24D22	Ø 2.4mm x 2200mm

18
Accessories
Valve
Cleaning Brushes



Double-ended valve cleaning brush
single-use, length 150mm
Ø brush head 5 and 11.5mm

EA-S1001PCD

Double-ended valve cleaning brush

unit = 10 pcs




Single-ended valve cleaning brush
single-use

EA-S1001 PC

Valve cleaning brush / single

unit = 10 pcs



Reusable valve cleaning brush

EA-R10PC

Valve cleaning brush, reusable

19
Accessories
Valve Caps



**Valve cap, single-use, compatible with
Olympus endoscope**

EA-SVC

valve cap

unit = 10 pcs



**Valve cap set, single-use, compatible with
Olympus und Fuji endoscope**

EA-SVC-S

valve cap set

Endoaccess GmbH
Feldriethe 1
D-30826 Garbsen
Phone: +49 5131-44 22 61 0
info@endoaccess.de
www.endoaccess.de